



# Recombinant YFV Nonstructural protein NS2A (Fragment) [His] (DAGA-1079)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Recombinant 6xHis tagged YFV Nonstructural protein NS2A (Fragment) (UniProtKB No. P89916) was expressed.
<b>Purity</b>	>85%
<b>Conjugate</b>	His
<b>Applications</b>	ELISA, Immunogen
<b>Recommended Usage</b>	This item requires custom production and lead time is between 4-8 weeks. We can custom produce based on your specifications.
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg, 5 mg
<b>Buffer</b>	PBS, pH7.4
<b>Preservative</b>	None
<b>Storage</b>	Store at 2-8 °C short term and -20°C long term. It is recommended to aliquot needed volumes and store at -20°C. Avoid repeated freeze/thaw cycles.

## BACKGROUND

**Keywords** Yellow fever virus; NS2A; Nonstructural protein

# GENE INFORMATION

UniProt ID	P89916
------------	--------

---