



Vitamin A [HRP] (DAGA-119H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Vitamin A (retinol) and its metabolites play diverse roles in physiology, ranging from incorporation into vision pigments to controlling transcription of a host of important genes.
Conjugate	HRP
Applications	ELISA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Retinal; Retinoic acid; Retinol;
----------	----------------------------------