



Vitamin E [KLH] (DAGA-128K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Vitamin E refers to a group of ten lipid-soluble compounds that include both tocopherols and tocotrienols. Of the many different forms of vitamin E, γ -tocopherol is the most common in the North American diet. γ -Tocopherol can be found in corn oil, soybean
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Vitamin E; Alpha-Tocopherol (vitamin E)
-----------------	---