



Androstenedione [KLH] (DAGA-133K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	In addition to its role as a precursor to testosterone in the body, androstenedione has androgenic properties in its own right, acting as a weak partial agonist of the androgen receptor. However, in the presence of full agonists like testosterone or dihydrotestosterone (DHT), due to its lower intrinsic activity in comparison, it has antagonistic properties, and can behave more like an antiandrogen like cyproterone.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Androstenedione; 4 androstenedione; 17 ketoestosterone; Androstanedione 3;
-----------------	--