



17b Boldenone [BSA] (DAGA-136B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Boldenone was synthesized in an attempt to create a long-acting injectable methandrostenolone (Dianabol), but in reality boldenone acts nothing like methandrostenolone. It has a very long half-life, and can show up on a steroid test for up to 1.5 years, due to the long undecylenate ester attached to the parent steroid. Trace amounts of the drug can easily be detected for months after discontinued use.
Conjugate	BSA
Applications	LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	17b Boldenone; 1 4-androstadiene-3-one-17b-ol; Equipoise; Ganabol; Equigan; Ultragan; Boldenone;
-----------------	--