



17b Boldenone [KLH] (DAGA-136K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Boldenone was synthesized in an attempt to create a long-acting injectable methandrostenolone (Dianabol), but in reality boldenone acts nothing like methandrostenolone. It has a very long half-life, and can show up on a steroid test for up to 1.5 years, due to the long undecylenate ester attached to the parent steroid. Trace amounts of the drug can easily be detected for months after discontinued use.

Conjugate KLH

Applications Immunogen

Format Liquid

Concentration Batch dependent - please inquire should you have specific requirements.

Size 1 mg

Buffer 0.01M pH7.4 PBS

Preservative None

Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords 17b Boldenone; 1 4-androstadiene-3-one-17b-ol; Equipoise; Ganabol; Equigan; Ultragan; Boldenone;