



B-Zeranol [HRP] (DAGA-138H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------------|--|
| Antigen Description | Zeranol is a synthetic non-steroidal anabolic agent, with oestrogenic activity. The main metabolites of zeranol in mammals are zeralanone and taleranol. Zeranol and its metabolites are excreted as free compounds and as glucuronide or sulphate tracer. Zeranol is also produced naturally from the mycotoxin zearalenone, which is found in seven strains of Fusaria that have been identified in animal feedstuffs. Therefore, any tissue residues of zeranol may be of natural origin. |
| Conjugate | HRP |
| Applications | ELISA |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Size | 1 mg |
| Buffer | 0.01M pH7.4 PBS |
| Preservative | None |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

BACKGROUND

Keywords beta Zeranol; B-Zeranol