



CLAD [BSA] (DAGA-140B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CLAD is a synthetic anabolic androgenic steroid.
Conjugate	BSA
Applications	LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	CLAD; 4 Androsten 4 Chloro 3 17 Dione; 4-Chlorotestosterone; Clostebol; Clostebolum;
----------	--