



Deoxy-Corticosterone [KLH] (DAGA-143K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	A mineralocorticoid produced in small quantities by the human adrenal cortex, having about 3 per cent of the sodium-retaining activity of aldosterone.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords cortexone; Deoxycorticosterone; deoxycortone; desoxycorticosterone; desoxycortone;