



Dihydrotestosterone [KLH] (DAGA-146K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Dihydrotestosterone is a biologically active metabolite of the hormone testosterone, formed primarily in the prostate gland, testes, hair follicles, and adrenal glands by the enzyme 5α-reductase by means of reducing the 4,5 double-bond.
----------------------------	--

Conjugate	KLH
------------------	-----

Applications	Immunogen
---------------------	-----------

Format	Liquid
---------------	--------

Concentration	Batch dependent - please inquire should you have specific requirements.
----------------------	---

Size	1 mg
-------------	------

Buffer	0.01M pH7.4 PBS
---------------	-----------------

Preservative	None
---------------------	------

Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
----------------	---

BACKGROUND

Keywords	5αDHT; DHT; INN;
-----------------	------------------