



Epitestosterone [BSA] (DAGA-147B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Testosterone is a steroid hormone from the androgen group. Testosterone is primarily secreted in the testes of males and the ovaries of females although small amounts are secreted by the adrenal glands.
Conjugate	BSA
Applications	LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	ISOTESTOSTERONE; DELTA-ANDROSTEN-17A-OL-3-ONE; CIS-TESTOSTERONE; EPITESTOSTERONE; 4-ANDROSTEN-17ALPHA-OL-3-ONE; 17a-testosterone; 17-hydroxy-androst-4-en-3-one; 17-epitestosterone; 17ALPHA-TESTOSTERONE
-----------------	---