



## Epitestosterone [KLH] (DAGA-147K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Testosterone is a steroid hormone from the androgen group. Testosterone is primarily secreted in the testes of males and the ovaries of females although small amounts are secreted by the adrenal glands.
<b>Conjugate</b>	KLH
<b>Applications</b>	Immunogen
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

<b>Keywords</b>	ISOTESTOSTERONE; DELTA-ANDROSTEN-17A-OL-3-ONE; CIS-TESTOSTERONE; EPITESTOSTERONE; 4-ANDROSTEN-17ALPHA-OL-3-ONE; 17a-testosterone; 17-hydroxy-androst-4-en-3-one; 17-epitestosterone; 17ALPHA-TESTOSTERONE
-----------------	---