



Ethinyl Estradiol 3 [KLH] (DAGA-150K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Ethinyl Estradiol is a derivative of 17b-estradiol (E2), the major endogenous estrogen in humans.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Ethinyl Estradiol 3; Ethinylestradiol; Ethynyl estradiol; Ethinyl estradiol; Estradiol;
----------	---