



Ethynylestradiol [BSA] (DAGA-152B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Ethinyl estradiol (EE), also sometimes written as ethinylestradiol, ethynyl estradiol, or ethinyl estradiol, is a derivative of 17 β -estradiol (E2), the major endogenous estrogen in humans.
Conjugate	BSA
Applications	LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords EE; Ethynylestradiol