



Flumethasone [KLH] (DAGA-153K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Flumethasone is a glucocorticoid. It is an anti-inflammatory and a corticosteroid for topical use. Besides, It is also used for the treatment of otitis externa and otomycosis.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Flumethasone; Flumetasone Pivalate; Locacorten; Locorten;
-----------------	---