



Gibberellic Acid [KLH] (DAGA-155K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Gibberellic acid (also called Gibberellin A3, GA, and GA3) is a hormone found in plants. Its chemical formula is C ₁₉ H ₂₂ O ₆ . When purified, it is a white to pale-yellow solid.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Gibberellic acid; Gibberellin A3; GA; GA3
----------	---