



# Hydroxyprogesterone [BSA] (DAGA-157B)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Antigen Description	Hydroxyprogesterone is a progestogen structurally related to progesterone that has been used for recurrent miscarriage and various menstrual disorders.
Conjugate	BSA
Applications	LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

Keywords	17 alpha hydroxyprogesterone caproate; 17 hydroxy progesterone; Progesterone 17OH; Progesterone8(17OH);
----------	---