



## Indole 3 Acetic Acid [KLH] (DAGA-159K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|                            |   |
|----------------------------|---|
| <b>Antigen Description</b> | Indole-3-acetic acid (IAA), is a member of the group of phytohormones called auxins and is generally considered to be the most important native auxin. It is synthesized by plants from the amino acid tryptophan. IAA also promotes growth of vascular tissue, |
| <b>Conjugate</b>           | KLH   |
| <b>Applications</b>        | Immunogen   |
| <b>Format</b>              | Liquid  |
| <b>Concentration</b>       | Batch dependent - please inquire should you have specific requirements.   |
| <b>Size</b>                | 1 mg  |
| <b>Buffer</b>              | 0.01M pH7.4 PBS   |
| <b>Preservative</b>        | None  |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.   |

### BACKGROUND

|                 |  |
|-----------------|--|
| <b>Keywords</b> | Indole 3 acetic acid; alpha-Indol-3-yl-acetic acid; alpha-indol-3-yl-aceticacid; beta-IAA; beta-Indole-3-acetic acid; beta-indole-3-aceticacid |
|-----------------|--|