



Isoxsuprine [KLH] (DAGA-161K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Isoxsuprine is most commonly used to treat hoof-related problems in the horse, most commonly for laminitis and navicular disease, as its effects as a vasodilator are thought to increase circulation within the hoof to help counteract the problems associate
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords Isoxsuprine; isoxsuprine hydrochloride; Vasodilan;