



Recombinant Chikungunya E2 Antigen (DAGA-162)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Chikungunya is an infection caused by the chikungunya virus. The disease features the sudden onset of fever two to four days after exposure. The fever usually lasts two to seven days, while accompanying joint pains typically last weeks or months but sometimes years. The mortality rate is a little less than 1 in 1,000; the elderly or those with underlying chronic medical problems are most likely to have severe complications.

Purity ≥95% as determined by SDS-PAGE

Conjugate TBD

Applications WB standard, antibody ELISA, immunogen, etc.

Molecular Weight 45.7 kDa

Concentration Batch dependent - please inquire should you have specific requirements.

Size 1 mg

Buffer 10mM PBS, pH7.4

Preservative None

BACKGROUND

Keywords CHIKV;Chikungunya Virus;Chikungunya E1;Chikungunya E2;CHIKV E1;CHIKV E2