



# Medroxyprogesterone Acetate [KLH] (DAGA-162K)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

**Antigen Description** Medroxyprogesterone acetate is a steroidal progestin, a synthetic variant of the human hormone progesterone. It is used as a contraceptive, in hormone replacement therapy and for the treatment of endometriosis as well as several other indications.

**Conjugate** KLH

**Applications** Immunogen

**Format** Liquid

**Concentration** Batch dependent - please inquire should you have specific requirements.

**Size** 1 mg

**Buffer** 0.01M pH7.4 PBS

**Preservative** None

**Storage** Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Keywords** Medroxyprogesterone Acetate; 17a-hydroxy-6a-methylprogesterone acetate; MPA;