



Melengestrol [BSA] (DAGA-164B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------|---|
| Antigen Description | Melengestrol is a progestogen-like compound. |
| Conjugate | BSA |
| Applications | LFIA |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Size | 1 mg |
| Buffer | 0.01M pH7.4 PBS |
| Preservative | None |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

BACKGROUND

| | |
|----------|-------------------------------------|
| Keywords | Melengestrol; Melengestrol acetate; |
|----------|-------------------------------------|