



Ritodrine [BSA] (DAGA-173B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Ritodrine, (4-[(1R,2S)-1-hydroxy-2-[2-(4-hydroxyphenyl)ethylamino]propyl]phenol), is an adrenergic beta-agonist used to control premature labour. It has a molecular weight of 287.35354 g/mol and molecular formula C₁₇H₂₁NO₃.

Conjugate BSA

Applications LFIA

Format Liquid

Concentration Batch dependent - please inquire should you have specific requirements.

Size 1 mg

Buffer 0.01M pH7.4 PBS

Preservative None

Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords 4 [(1R; 2S) 1 hydroxy3[2 (4 hydroxyphenyl)ethylamino]propyl]phenol; Ritodrina; Ritodrinium; Yutopar;