



Salicylic Acid [KLH] (DAGA-174K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description Salicylic acid is a monohydroxybenzoic acid, a type of phenolic acid and a beta hydroxy acid. This colorless crystalline organic acid is widely used in organic synthesis and functions as a plant hormone. It is derived from the metabolism of salicin. In addition to being a compound that is chemically similar to but not identical to the active component of aspirin (acetylsalicylic acid), it is probably best known for its use in anti-acne treatments.

Conjugate KLH

Applications Immunogen

Format Liquid

Concentration Batch dependent - please inquire should you have specific requirements.

Size 1 mg

Buffer 0.01M pH7.4 PBS

Preservative None

Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords Salicylic Acid; C₆H₄(OH)COOH; Aspirin;