



## Salmeterol [BSA] (DAGA-175B)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

**Antigen Description** Salmeterol is a long-acting beta2-adrenergic receptor agonist drug that is currently prescribed for the treatment of asthma and chronic obstructive pulmonary disease(COPD). It is currently available as a dry powder inhaler that releases a powdered form of the drug. Before 2008, it was also available as a metered-dose inhaler (MDI).

**Conjugate** BSA

**Applications** LFIA

**Format** Liquid

**Concentration** Batch dependent - please inquire should you have specific requirements.

**Size** 1 mg

**Buffer** 0.01M pH7.4 PBS

**Preservative** None

**Storage** Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

**Keywords** beta2-adrenergic receptor agonist drug;