



## Stanozolol [BSA] (DAGA-176B)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|                            |   |
|----------------------------|---|
| <b>Antigen Description</b> | Stanozolol has been used in both animal and human patients for a number of conditions. In humans, it has been demonstrated to be successful in treating anaemia and hereditary angioedema. Veterinarians may prescribe the drug to improve muscle growth, red blood cell production, increase bone density and stimulate the appetite of debilitated or weakened animals. |
| <b>Conjugate</b>           | BSA   |
| <b>Applications</b>        | LFIA  |
| <b>Format</b>              | Liquid  |
| <b>Concentration</b>       | Batch dependent - please inquire should you have specific requirements.   |
| <b>Size</b>                | 1 mg  |
| <b>Buffer</b>              | 0.01M pH7.4 PBS   |
| <b>Preservative</b>        | None  |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.   |

### BACKGROUND

|                 |                                       |
|-----------------|---------------------------------------|
| <b>Keywords</b> | Stanozolol; Winstrol; Winstrol Depot; |
|-----------------|---------------------------------------|