



Trenbolone [KLH] (DAGA-178K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Trenbolone is a steroid used by veterinarians on livestock to increase muscle growth and appetite. To increase its effective half-life, trenbolone is not used in an unrefined form, but is rather administered as ester derivatives such as trenbolone acetate.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords Trenbolone; 17 β -Hydroxyestra-4; 9; 11-trien-3-one