



17a Trenbolone [KLH] (DAGA-179K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Trenbolone is a steroid used by veterinarians on livestock to increase muscle growth and appetite. To increase its effective half-life, trenbolone is not used in an unrefined form, but is rather administered as ester derivatives such as trenbolone acetate, trenbolone enanthate or trenbolone cyclohexylmethylcarbonate (Parabolan). Plasma lipases then cleave the ester group in the bloodstream leaving free trenbolone.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords 17a Trenbolone; 17b-Hydroxyestra-4 9 11-trien-3-one; Trenbolone;