



Triamcinolone [HRP] (DAGA-183H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Triamcinolone Acetonide is a more potent type of triamcinolone (synthetic corticosteroid), being about 8 times as effective as prednisone. It is used to treat medical conditions, such as eczema, psoriasis, arthritis, allergies, ulcerative colitis, lupus, sympathetic ophthalmia, temporal arteritis, uveitis, and ocular inflammation, visualization during vitrectomy and the prevention of asthma attacks. Triamcinolone Acetonide has an elimination half-life of 1.5-2.5 hours in Humans.
Conjugate	HRP
Applications	ELISA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Triamcinolone Acetonide; Kenalog; Polcortolon; Aristoderm; Aristogel; Triaceton; Tricinelon; Azmacort; kenalone; Nasacort; Solodelf;
-----------------	--