



## Triiodothyronine [KLH] (DAGA-184K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Triiodothyronine, C <sub>15</sub> H <sub>12</sub> I <sub>3</sub> NO <sub>4</sub> , also known as T <sub>3</sub> , is a thyroid hormone. It affects almost every physiological process in the body, including growth and development, metabolism, body temperature, and heart rate.
<b>Conjugate</b>	KLH
<b>Applications</b>	Immunogen
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

**Keywords** triiodothyronine; reverse; TRIIODOTHYRONINE; REVERSE T<sub>3</sub>; T<sub>3</sub>; LIOTHYRONIN