



Erythromycin [KLH] (DAGA-189K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Erythromycin is amacrolide antibiotic that has an antimicrobial spectrum similar to or slightly wider than that of penicillin, and is often used for people who have an allergy to penicillins.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Erythromycin; Actinobacteria; Actinobacteria; Actinobacteridae; Actinomycetales; Pseudonocardia; Saccharopolyspora; S. erythraea; Saccharopolyspora erythraea
-----------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------