



Flumequine [BSA] (DAGA-193B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Flumequine is a first generation drug of the fluoroquinolone group of antibiotics. It has general use in the treatment of enteric infections in domestic species.
Conjugate	BSA
Applications	LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	9-Fluoro-6; 7-dihydro-5-methyl-1-oxo-1H; 5H-benzo[ij]-quinolizine-2-carboxylic acid; Apurone;
-----------------	---