



## Ormetoprim [OVA] (DAGA-1970)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	10mM PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

### BACKGROUND

Introduction	Ormetoprim is a close analog of diaveridine and a folic acid antagonist used in potentiated sulfonamide mixtures used as coccidiostats.
--------------	---