



## Ronidazole [BSA] (DAGA-198B)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Ronidazole is an antiprotozoal agent used in veterinary medicine for the treatment of histomoniasis and swine dysentery. It may also have use for the treatment of Tritrichomonas foetus infection in cats.
<b>Conjugate</b>	BSA
<b>Applications</b>	LFIA
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

<b>Keywords</b>	Ronidazole
-----------------	------------