



Sulfachloropyrazine [BSA] (DAGA-201B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The sulphonamides are distributed widely throughout all organs and tissues; and are metabolized primarily by the liver, producing inactive acetylated and glucuronide forms, which are excreted in urine.
Conjugate	BSA
Applications	LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords Sulphachloropyrazine; Sulfachloropyrazine Sodium;