



# Sulfachloropyrazine [KLH] (DAGA-201K)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	The sulphonamides are distributed widely throughout all organs and tissues; and are metabolized primarily by the liver, producing inactive acetylated and glucuronide forms, which are excreted in urine.
<b>Conjugate</b>	KLH
<b>Applications</b>	Immunogen
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Keywords** Sulphachloropyrazine; Sulfachloropyrazine Sodium;