



## Sulfadoxine [BSA] (DAGA-202B)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Sulfadoxine (also spelled sulphadoxine) is an ultra-long-lasting sulfonamide often used in combination with pyrimethamine to treat or prevent malaria. It is also used, usually in combination with other drugs, to treat or prevent various infections in livestock.
<b>Conjugate</b>	BSA
<b>Applications</b>	LFIA
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

<b>Keywords</b>	sulphadoxine; sulphormethoxine; wr4873; SULFADOXIN; SULFADOXINE; sulformethoxine; sulforthomidine; 6-(4-aminobenzenesulfonamido)-4; 5-dimethoxypyrimidine; 4-sulfanilamido-5; 6-dimethoxypyrimidine
-----------------	---