



Sulfaguanidine [KLH] (DAGA-203K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Sulfaguanidine is a sulfonamide.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	p-Aminobenzenesulfoguanidide; Resulfon; RP 2275; rp2275; Ruocid; SG; S-Guanidine; Shigatox; Suganyl; Sulfaguamide; Sulfaguanidin; Sulfaguanil; Sulfaguine; Sulfanilamide; N(sup1)-amidino-; Sulfanilamide; N1-(diaminomethylene)-; Sulfanilamide; N1-amidino-; Sulfanilguanidine; sulfanilyl-guanidin; Sulfentidine; Sulfoguanidine; Sulfoguanil; Sulfoguanyl; Sulfoguenil; Sulfoquanidine;
----------	---