



## Sulfapyridine [BSA] (DAGA-208B)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

**Antigen Description** Sulphapyridine, 4-amino-N-pyridin-2-yl-benzenesulfonamide, is an antibacterial agent used to treat certain skin diseases. It is a competitive inhibitor of bacterial para-aminobenzoic acid (PABA), a substrate of the enzyme dihydropteroate synthetase.

**Conjugate** BSA

**Applications** LFIA

**Format** Liquid

**Concentration** Batch dependent - please inquire should you have specific requirements.

**Size** 1 mg

**Buffer** 0.01M pH7.4 PBS

**Preservative** None

**Storage** Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### BACKGROUND

**Keywords** 2 Sulfapyridine; 4 amino N pyridin3yl benzenesulfonamide; Coccoclase; Eubasinum; Haptocil; Piridazol; Relbapiridina; Septipulmon; sulfapyridine; Thioseptal;