



## Sulfapyridine [KLH] (DAGA-208K)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Sulphapyridine, 4-amino-N-pyridin-2-yl-benzenesulfonamide, is an antibacterial agent used to treat certain skin diseases. It is a competitive inhibitor of bacterial para-aminobenzoic acid (PABA), a substrate of the enzyme dihydropteroate synthetase.
----------------------------	---

<b>Conjugate</b>	KLH
------------------	-----

<b>Applications</b>	Immunogen
---------------------	-----------

<b>Format</b>	Liquid
---------------	--------

<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
----------------------	---

<b>Size</b>	1 mg
-------------	------

<b>Buffer</b>	0.01M pH7.4 PBS
---------------	-----------------

<b>Preservative</b>	None
---------------------	------

<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
----------------	---

### BACKGROUND

<b>Keywords</b>	2 Sulfapyridine; 4 amino N pyridin3yl benzenesulfonamide; Coccoclase; Eubasinum; Haptocil; Piridazol; Relbapiridina; Septipulmon; sulfapyridine; Thioseptal;
-----------------	--