



Tobramycin [KLH] (DAGA-213K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Tobramycin is an aminoglycoside antibiotic used to treat various types of bacterial infections, particularly Gram-negative infections.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Tobramycin; Anti – Tobramycin; 1-epitobramycin; 3'-deoxykanamycinb; d-6-tyrideoxy-alpha-d-ribohexopyranosyl-(1-6))-2-deoxy; distobram gernebcin lpha-d-ribohexopyranosyl-(1-4))-2-deoxy-d-streptamine; nebcin; nebicina; nebramycin6; nebramycinvi
-----------------	--