



Acepromazine [KLH] (DAGA-217K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Acepromazine is a phenothiazine derivative antipsychotic drug. It was first used in humans in the 1950s,[1] but is now little used in humans (the closely related analogue, chlorpromazine, is still used as an antipsychotic in humans). Acepromazine is frequ
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Acepromazine
-----------------	--------------