



Acetaminophen [KLH] (DAGA-218K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ACETAMINOPHEN is a widely used over-the-counter analgesic (pain reliever) and antipyretic (fever reducer). It is commonly used for the relief of headaches and other minor aches and pains and is a major ingredient in numerous cold and flu remedies. In combination with opioid analgesics, paracetamol can also be used in the management of more severe pain such as post-surgical pain and providing palliative care in advanced cancer patients. The onset of analgesia is approximately 11 minutes after oral administration of paracetamol, and its half-life is 1–4 hours. Though acetaminophen is used to treat inflammatory pain, it is not generally classified as an NSAID because it exhibits only weak anti-inflammatory activity.
----------------------------	--

Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	ACETAMINOPHEN; N-(4-hydroxyphenyl)ethanamide; N-(4-hydroxyphenyl)acetamide;
-----------------	---