



TSH protein (>95% pure) (DAGA-225)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95 % as determined by SDS-PAGE
Conjugate	N/A
Applications	Standard
Bio-activity	≥ 6.2 IU/mg
Format	Lyophilized
Size	10 µg, 50 µg, 100 µg, 1 mg
Buffer	10 mM Sodium Phosphate, 150 mM Sodium Chloride, 1 mg/ml BSA, pH 7.4
Preservative	0.1% Sodium Azide (optional)
Storage	Store at 2°C to 8°C short term (up to 2 weeks) or frozen below -20°C

BACKGROUND

Introduction	Thyroid-stimulating hormone (also known as TSH or thyrotropin) is a peptide hormone synthesized and secreted by thyrotrope cells in the anterior pituitary gland, which regulates the endocrine function of the thyroid gland.
Keywords	Thyroid stimulating hormone; TSH