



# Amphetamine [KLH] (DAGA-225K)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Amphetamine (USAN) or amphetamine (INN) is a psychostimulant drug of the phenethylamine class which produces increased wakefulness and focus in association with decreased fatigue and appetite. Brand names of medications that contain, or metabolize into, amphetamine include Adderall, Dexedrine, Dextrostat, Desoxyn, ProCentra, and Vyvanse, as well as Benzedrine in the past.
<b>Conjugate</b>	KLH
<b>Applications</b>	Immunogen
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	1 mg
<b>Buffer</b>	0.01M pH7.4 PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

<b>Keywords</b>	Propylamine; alpha-methyl-gamma-phenyl-; Propylamine; .alpha.-methyl-.gamma.-phenyl-; (+/-)-AMPHETAMINE; DL-AMPHETAMINE; AURORA KA-7110; 1-Phenyl-3-aminobutane; 3-Phenyl-1-methylpropylamine; 4-Phenyl-2-aminobutane; Amphetamine; AMP; (n)-3-amino-1-phenylbutane; (n)-4-phenyl-2-butylamine; METHYL-3-PHENYLPROPYLAMINE;
-----------------	---