



Azaperone [KLH] (DAGA-227K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Azaperone is a sedative and antiemetic, which is used mainly as a tranquilizer in veterinary medicine to reduce aggression in farmed pigs. More rarely it may be used in humans as an antipsychotic drug, but this is uncommon as it may cause hypotension or respiratory depression.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Stresnil; Fluoperidol; Azapero; Azaperone
-----------------	---