



Luteinizing Hormone (partially purified) (DAGA-227P)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	~40%
Conjugate	N/A
Applications	N/A
Format	Lyophilized
Size	10IU
Buffer	Lyophilized from 20mM PB, pH 7.4
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction Luteinizing hormone is a hormone produced by the anterior pituitary gland. In females, an acute rise of LH called the LH surge triggers ovulation and development of the corpus luteum. In males, where LH had also been called interstitial cell-stimulating hormone (ICSH), it stimulates Leydig cell production of testosterone.

Keywords Luteinizing hormone; LH