



Carboxymethyl-lysine [BSA] (DAGA-236B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	N(7)-Carboxymethyllysine (CML), also known as N(epsilon)-(carboxymethyl)lysine, is an advanced glycation endproduct (AGE). CML has been the most used marker for AGEs in food analysis.
Conjugate	BSA
Applications	LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Carboxymethyl-lysine
----------	----------------------