



Clonidine [KLH] (DAGA-238K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Clonidine is an alpha 2-adrenoceptor agonist. Clonidine is extensively used as an adjunct to spinal morphine and is suggested to be particularly effective for neuropathic pain; its clinical use is limited by side effects such as sedation and hypotension.
Conjugate	KLH
Applications	Immunogen
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Clonidine; alpha 2-adrenoceptor agonist
-----------------	---