



Clonidine [OVA] (DAGA-238O)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LFIA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	10mM PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	Clonidine is an alpha 2-adrenoceptor agonist. Clonidine is extensively used as an adjunct to spinal morphine and is suggested to be particularly effective for neuropathic pain; its clinical use is limited by side effects such as sedation and hypotension.
--------------	--